

S S Technomed (P) Limited
(AN ISO 9001:2008 & 13485-2003 CERTIFIED COMPANY)



India's Leading Manufacturer of Pediatric Equipments



INFANT RADIANT WARMER
(MODEL : TIANA™ - NW)

INFANT RADIANT WARMER (MODEL : TIANA™ - NW)

Advance micro processor controlled system with self check on the electronic circuitry Skin / Manual control. Digital Display of Baby Temperature / Set Temperature. Manual / Servo two function control modes provide the clinician more flexibility while optimizing the thermal environment. The Spacious baby bed incorporated easily drop down sides for maximum access to the baby. The Control Panel is also the heart for a completely integrated alarm system with both audible and visual Indicators in the event of

- Probe Fail
- Skin High
- Skin Low
- Power Fail

Manual Mode

Facility of converting the warmer from servo mode to manual mode for non-stop working. The heat in the manual mode is locked at 50% to give more protection for the baby.

Heater Assembly

The Heater assembly located above the center of the support structure consists of radiant heater and reflector with inbuilt halogen examination lamp for assessment of Infant color. The heater can be swivelled to either direction.

Baby Bed

The attached spacious baby bed incorporates easy drop down glasses for maximum access to the patient.

Castors

The system is fitted with 2" high quality castors having brakes in front 2 castors. Antistatic castors can be fitted as an option for further safety in lieu of the standard castors.

Dimensions(Approx.)

Technical Specification

Height- 1920mm

Length- 890mm

Width- 860mm

Electrical Specifications

Power Supply 220V/50Hz

Heater Power 650 Watts

Fuse 5 Amp

OPTIONAL

2 Trays below the baby bed for keeping consumables.

2 Drawers fitted below the baby bed.

Technical Specifications are subject to change without prior notice. As manufacturers we reserve the right to change / modify any / all the specifications due to philosophy of continuous upgradation.

S S Technomed (P) Limited

(AN ISO 9001:2008 & 13485-2003 CERTIFIED COMPANY)

Works: A-128, Sector - A4, Tronica City, UPSIDC Industrial Area, Loni, Ghaziabad - 201102 Uttar Pradesh (INDIA)

Phone No: +91 - 120 - 2696390, Fax No: +91 - 120 - 2696391

E - mail: sstechnomed2@gmail.com , bhavna@sstechnomed.com

Website : www.sstechnomed.com

